

# → STARLIGHT UNO. COLORFUL AND EVEN EASIEST HANDLING.

New LED curing light, new colours!  
One touch handling, outstanding  
curing performance, exchangeable  
battery ... persuasive features,  
persuasive price!

## → STARLIGHT UNO

1 handpiece, 1 charging station

Color

- white
- grey
- turquoise
- yellow
- light blue
- violet

Ref. no.

- 05100092
- 05100093
- 05100094
- 05100095
- 05100096
- 05100097





available from June 2017

## → TECHNICAL DATA

**sizes (l x Ø in mm):**

190 x 21

**weight:**

77 g

**light source:**

high power LED

**wavelength:**

dominant wavelength 440 – 465 nm

**light intensity:**

> 1.500 mW/cm<sup>2</sup>

**average life span LED:**

1.800.000 cycles of 20 sec

**battery:**

rechargeable Lithium-ion battery  
(no memory-effect)

**capacity of the accumulator:**

240 cycles of 10 sec

**charging time:**

4 hours (240 min)

**warranty:**

3 years, battery 1 year

## → STARLIGHT UNO ACCESSORIES

Product	Ref. no.
<b>battery starlight uno</b>	03020222
<b>protective sleeves</b> (500 pieces)	03230015
<b>optical protection</b> (package with 2 pieces)	02900130