

Next generation of top dental innovations.

Easy-to-use X-ray with the highest image resolution. Designed for comfort. Large active area for maximum coverage. Try it.



Designed and manufactured in the world's leading facility of CCD and CMOS technology for guaranteed highest-grade imaging.

Dentalmind Digital X-ray II.

Size 1 and 2 sensor includes the following:

- 2-year limited warranty
- Sensor Cable Length: up to 20 meters
- Connects directly to your computer via USB 2.0
- Repairable USB cable for longevity of sensors
- Fiber optic plate sensors for low-noise ratio
- Fully waterproof with sealed cable and housing design

- Largest active area for maximum coverage
- High resolution, the image is exceptionally sharp with great clarity
- Sensor Technology: Next Generation CMOS Fiber optic Sensors
- Image Pixel Size: 19x19
- Maximum Gray Levels: 4096 (12-bit)
- CE marked in class IIa CE 0088

Size 1 sensor specifications

- Image Size - 1580 x 1050 pixels
- External Dimensions: 24.35 x 36.73 mm
- Active Area - 19.95 x 30.02 mm

Size 2 sensor specifications

- Image Size: 1896 x 1368 Pixels
- External Dimensions: 30.49 x 42.80 mm
- Active Area: 25.99 x 36.02 mm

www.dentalmind.com

dentalmind®

DIGITAL X-RAY II