

# Our intraoral camera is the high tech definition of clarity.

High end lens, hand control focus, and compatibility with almost all dental imaging software. Award winning.



Digital precision. Excellent image quality. Optics with focal range of 1 mm-infinity for macro to full facial shots. Easy handling.

## Dentalmind Intraoral Camera.

- Precision optics with focal range of 1 mm-infinity for macro to full facial shots
- Precision controlled sliding focus
- 8 white LED lights.
- Light ON/OFF capability
- Ergonomic hand-piece that weighs only 2.4 oz
- Silicon capture button for image stability
- Secure Monster® USB cable connection
- Integrates seamlessly with almost all dental imaging software.
- Imager: Sony Super HAD Color 1/4"
- Resolution: 720 X 480 Pixels
- Frame Rate: 30 Fps Max
- Software Interface: WDM
- Data Format: YUV 4:2:2. ITU-R 656
- Field of View: 60 Degrees
- Focus Type: Variable Focus
- Signal to Noise Ratio: 51 Db
- Video Output: USB 2.0 Digital
- Power Requirements: 5v Dc, Powered By USB Port
- Dimensions: 8.3(L) X 0.85(W) X 1.15(H) Inches
- 2-year Manufacturer Warranty, including parts and labor
- CE mark MDD 93/42 EEC

[www.dentalmind.com](http://www.dentalmind.com)

**dentalmind®**  
INTRAORAL CAMERA